

CERTIFICATE OF ANALYSIS


Reference Substance

Date

Cephalexin Impurity J

Batch No.:

IUPAC Name: (6R,7R)-7-((R)-2-amino-2-phenylacetamido)-3-methyl-8-oxo-5-thia-1-azabicyclo[4.2.0]oct-2-ene-2-carboxylic acid 5-oxide



Molecular Formula: C₁₆H₁₇N₃O₅S

Molecular Weight: 363.39

Sr.#	Test	Results	Data ref.
1	Appearance	Off white solid	-----
2	Purity by HPLC	98.61%	Chromatogram attached
3	MASS	Conform	Report attached
4	¹ H NMR	Conform	Report attached
5	IR	Conform	Report attached

Storage condition: Store at 2-8°C and protect from moisture and direct sunlight.

Shipment date:

This certificate of analysis is valid for two years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by