

CERTIFICATE OF ANALYSIS


Reference Substance

Date

Azilsartan K Medoxomil Impurity-5 (Bis Impurity)

Batch No.:

IUPAC Name: (5-methyl-2-oxo-1,3-dioxol-4-yl)methyl 2-ethoxy-1-((2'-(4-((5-methyl-2-oxo-1,3-dioxol-4-yl)methyl)-5-oxo-4,5-dihydro-1,2,4-oxadiazol-3-yl)-[1,1'-biphenyl]-4-yl)methyl)-1H-benzo[d]imidazole-7-carboxylate


 Molecular Formula: C₃₅H₂₈N₄O₁₁

Molecular Weight: 680.62

Sr.#	Test	Results	Data ref.
1	Appearance	Off white solid	-----
2	Purity by HPLC	93.29%	Chromatogram attached
3	MASS	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Storage condition: Store at 2-8°C and protect from moisture and direct sunlight.

Shipment date:

This certificate of analysis is valid for two years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:
Retest Date:
Document data reference:
Prepared by
Reviewed by
Approved by