

CERTIFICATE OF ANALYSIS


Reference Substance

Date

Ascorbic Acid

Batch No.

IUPAC Name: (R)-5-((S)-1, 2-dihydroxyethyl)-3,4-dihydroxyfuran-2(5H)-one

Molecular Formula: C₆H₈O₆

Molecular weight: 176.12 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off white Solid	-----
2	Purity by HPLC	96.53%	Chromatogram attached
3	¹ HNMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store at 2-8°C, protect from moisture.**Shipment date:**

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest date:****Document data reference:***Prepared by**Reviewed by**Approved by*