

## CERTIFICATE OF ANALYSIS


Reference Substance

Date

Butyric Acid

**Batch No.**

IUPAC Name: (1-Propanecarboxylic acid)



Molecular Formula:C4H8O2

Molecular weight:88.11 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Colorless Liquid	-----
2	Purity by GC	99.9%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform	Report attached

**Storage condition:** Store at Refrigerator, Protect from direct sunlight.**Shipment date:**

This certificate of analysis is valid for one year from date of shipment, provided the substance is stored under suitable conditions.

**Document data reference:***Prepared by**Reviewed by**Approved by*