

CERTIFICATE OF ANALYSIS


Reference Substance

Date:

Fexofenadine Impurity C

Batch No.:

IUPAC Name: (1RS)-4-[4-(hydroxydiphenylmethyl) piperidin-1-yl]-1-[4-(1-methylethyl)phenyl]butan-1-ol.


 Molecular Formula: $C_{31}H_{39}NO_2$

Molecular Weight: 457.65

Sr.#	Test	Results	Data ref.
1	Appearance	Light Brown solid powder	-----
2	Purity by HPLC	98.39%	Chromatogram attached
3	MASS	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Storage condition: Store at 2-8 degree temperature and protect from direct sunlight.

Shipment date:

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by