


CERTIFICATE OF ANALYSIS

Reference Substance

Date:

Dexamethasone- Δ 17,20 21-Aldehyde**Batch No.****IUPAC Name:** Name: 9-Fluoro-11 β ,20-dihydroxy-16 α -methyl-3-oxopregna-1,4,17(20)-trien-21-alMolecular Formula: C₂₂H₂₇FO₄

Molecular Weight: 374.46g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	White solid.	-----
2	Purity by HPLC	99.32%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached

Storage condition: Store at 2-8 degree, protect from direct sunlight and moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*