


## CERTIFICATE OF ANALYSIS

Reference Substance

Date:

**Desoxymetasone(17-Desoxy of Dexamethasone)****Batch No.****IUPAC Name:** Name: (11 $\beta$ ,16 $\alpha$ )-9-Fluoro-11,21-dihydroxy-16-methyl-pregna-1,4-diene-3,20-dione; Dexamethasone Impurity FMolecular Formula: C<sub>22</sub>H<sub>29</sub>FO<sub>4</sub>

Molecular Weight: 376.46g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	White solid.	-----
2	Purity by HPLC	98.25%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached

**Storage condition:** Store at 2-8 degree, protect from direct sunlight and moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

**Manufacturing date:****Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*