


CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard

Date:

9-DESFLUORO-11-KETO OF DEXAMETHASONE**Catalogue no.: VL273005**

Batch No:

IUPAC Name: (16 α)-17,21-Dihydroxy-16-methylpregna-1,4-diene-3,11,20-trioneMolecular Formula: C₂₂H₂₈O₅

Molecular weight: 372.45 g/mol

Sr. NO.	Test	Results	Data ref.
1	Appearance	White Solid	-----
2	Purity by HPLC	96.55%	Chromatogram attached
3	Mass	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest date:****Document data reference:**

Prepared by

Reviewed by

Approved by