


## CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard

Date:

**Dexamethasone 9,11-  $\alpha$  Epoxide****Catalogue no.: VL273006**

Batch No:

IUPAC Name: (11 $\alpha$ ,16 $\alpha$ )-9,11-Epoxy-17,21-dihydroxy-16-methylpregna-1,4-diene-3,20-dioneMolecular Formula: C<sub>22</sub>H<sub>28</sub>O<sub>5</sub>

Molecular weight: 372.45g/mol

Sr. NO.	Test	Results	Data ref.
1	Appearance	White Solid	-----
2	Purity by HPLC	99.46%	Chromatogram attached
3	Mass	Conform	Report attached
4	<sup>1</sup> H NMR	Conform	Report attached

**Material shipping condition:** At normal temperature.**Long term Storage:** Store at 2-8°C, Protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing provided the substance is stored under suitable conditions.

**Manufacturing date:****Retest date:****Document data reference:***Prepared by**Reviewed by**Approved by*