


CERTIFICATE OF ANALYSIS

Reference Substance

Date:

Bromfenac Sodium impurity- A**Batch No.:****IUPAC Name:** Sodium (2-amino-3-benzoylphenyl)acetate

Molecular Weight: 277.25

Molecular Formula: $C_{15}H_{12}NNaO_3$

Sr.#	Test	Results	Data ref.
1	Appearance	Yellow colored powder	-----
2	Purity by HPLC	98.31%	Chromatogram attached
3	Mass	Conform	Report attached
4	1H NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from moisture and direct sunlight.**Shipment date:**

This certificate of analysis is valid for two years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*