


CERTIFICATE OF ANALYSIS

Reference Substance

Date:

Bromfenac Sodium impurity- D**Batch No.:****IUPAC Name:** 7-benzoyl-1,3-dihydro-2H-indol-2-one

Molecular Weight: 237.25

Molecular Formula: C₁₅H₁₁NO₂

Sr.#	Test	Results	Data ref.
1	Appearance	Slightly Pink colored powder	-----
2	Purity by HPLC	99.40%	Chromatogram
3	Mass	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from moisture and direct sunlight.**Shipment date:**

This certificate of analysis is valid for two years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*