

CERTIFICATE OF ANALYSIS


Reference Substance

Date

Travoprost Related Compound A

Batch No.

IUPAC Name: (5Z)-7-[(1R, 2R, 3R,5S)-3,5-dihydroxy-2-[(1E,3R)-3-hydroxy-4-[3-(trifluoromethyl)phenoxy]-1-butenyl]cyclopentyl]-5-heptenoic Acid



Molecular Formula: C₂₃H₂₉F₃O₆

Molecular weight: 458.47 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Colorless Semisolid	-----
2	Purity by HPLC	96.71%	Chromatogram attached
3	MASS	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Storage condition: Store at 2-8 °C, protect from moisture and direct sunlight.

Shipment date:

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by