

CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD


Reference Standard
Clopidogrel (R) & (S) Isomer Mixture

Date:

Catalogue no.: VL31005

Batch No.:

IUPAC Name: (RS)-methyl 2-(2-chlorophenyl)-2-(6,7-dihydrothieno[3,2-c]pyridin-5(4H)-yl)acetate



Molecular Formula: C₁₆H₁₆ClNO₂S

Molecular Weight: 321.82 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	White Solid	-----
2	Purity by HPLC	97.26%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached

Storage condition: Store at 2-8 degree, protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by