

CERTIFICATE OF ANALYSIS


Reference Substance

Date

Cetirizine Impurity A/ Chlorobenzhydryl piperazine (CBHP)

Batch No.:

IUPAC Name: 1-((4-chlorophenyl)(phenyl)methyl)piperazine

Molecular Formula: C₁₇H₁₉ClN₂

Molecular weight: 286.80g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off White solid	-----
2	Purity by HPLC	99.95%	Chromatogram attached
3	Mass	Conform	Report attached
4	NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight and moisture.**Shipment date:**

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest date:****Document data reference:***Prepared by**Reviewed by**Approved by*