


CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard

Date:

Budesonide EP impurity L/Budesonide related compound L/11-Ketobudesonide**Catalogue No.:** VL34007

Batch No.

IUPAC Name: (16 α)-16,17-[Butylidenebis (oxy)]-21-hydroxypregna-1,4-diene-3,11,20-trioneMolecular Formula: C₂₅H₃₂O₆

Molecular weight: 428.52 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off White Solid Powder	-----
2	Purity by HPLC	>95.00%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight and Moisture.

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest date:****Document data reference:***Prepared by**Reviewed by**Approved by*