

CERTIFICATE OF ANALYSIS


Reference Substance

Date

Zolmitriptan Amine

Batch No.

IUPAC Name: (4S)-4-{{[3-(2-aminoethyl)-1H-indol-5-yl]methyl}}-1,3-oxazolidin-2-one

Molecular Formula: C₁₄H₁₇N₃O₂

Molecular weight: 259.30 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off-white colored solid	-----
2	Purity by HPLC	97.34%	Chromatogram attached
3	Mass	Conform	Report attached
4	NMR	Conform	Report attached

Storage condition: Store at 2-8°C, Protect from direct sunlight and moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest date:****Document data reference:***Prepared by**Reviewed by**Approved by*