

**CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD**


Reference Standard

Date

**Zolmitriptan R-Isomer****Catalog No.: VL350005**

Batch No.

IUPAC Name: (R)-4-[[3-[2-(Dimethylamino)ethyl]indol-5-yl]methyl]-2-oxazolidinone

Molecular Formula: C<sub>16</sub>H<sub>21</sub>N<sub>3</sub>O<sub>2</sub>

Molecular weight: 287.36 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off White Powder	-----
2	Purity by HPLC	98.70 %	Chromatogram attached
3	Mass	Conform	Report attached
4	NMR	Conform	Report attached
5	FTIR	Conform	Report attached

**Storage condition:** Store at 2-8°C, Protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

**Date of manufacturing:****Retest date:****Document data reference:***Prepared by**Reviewed by**Approved by*