

CERTIFICATE OF ANALYSIS


Reference Substance

Date:

Fesoterodine Impurity-B (Fesoterodine Impurity)

Batch No.:

IUPAC Name: - (R)-2-(3-(diisopropylamino)-1-phenylpropyl)-4-((isobutyryloxy)methyl)phenyl isobutyrate

Molecular Formula: C₃₀H₄₃NO₄

Molecular Weight: 481.68

Sr.#	Test	Results	Data ref.
1	Appearance	Brown liquid	-----
2	Purity by HPLC	96.21%	Chromatogram attached
3	¹ HNMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store at 2-8 degree, protect from direct sunlight & Moisture

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*