

## CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD


Reference Standard

Date:

**Fesoterodine Diol Dimer Impurity**

Batch No.

IUPAC Name: 2-((R)-3-(diisopropylamino)-1-phenylpropyl)-4-(((3-((S)-3-(diisopropylamino)-1-phenylpropyl)-4-hydroxybenzyl) oxy) methyl)phenol



Molecular Weight: 664.98

Molecular Formula: C<sub>44</sub>H<sub>60</sub>N<sub>2</sub>O<sub>3</sub>

Sr.#	Test	Results	Data ref.
1	Appearance	Yellow Sticky Liquid	-----
2	Purity by HPLC	91.06%	Chromatogram attached
3	Mass	Conform	Report attached
4	NMR	Conform	Report attached

**Storage condition:** Store at 2-8°C, protect from direct sunlight and moisture.

This certificate of analysis is valid for two years from date of manufacturing provided the substance is stored under suitable conditions.

**Date of manufacturing:**

**Retest Date:**

**Document data reference:**

*Prepared by*

*Reviewed by*

*Approved by*