


CERTIFICATE OF ANALYSIS

Reference Substance

Date:

Desloratadine N-Hydroxypiperidine / N-Hydroxy Derivative of Desloratidine**Batch No.:****IUPAC Name: Name:** 8-Chloro-6,11-dihydro-11-(1-hydroxy-4-piperidinylidene)-5H-benzo[5,6]cyclohepta[1,2-b]pyridineMolecular Formula: C₁₉H₁₉ClN₂O

Molecular Weight: 326.82 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Light Brown Solid	-----
2	Purity by HPLC	98.47%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached

Storage condition: Store at 2-8 degree, protect from direct sunlight and moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*