

## CERTIFICATE OF ANALYSIS


Reference Substance

Date

**Dipyridamole Impurity -A**

Batch No.

IUPAC Name: 2,2'-[[4,6,8-tri(piperidin-1-yl)pyrimido[5,4-d]pyrimidin-2-yl]nitrido]diethanol

Molecular Formula: C<sub>25</sub>H<sub>40</sub>N<sub>8</sub>O<sub>2</sub>

Molecular weight: 484.64 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Yellow colour solid	-----
2	Purity by HPLC	95.36%	Chromatogram attached
3	Mass	Conform	Report attached
4	<sup>1</sup> H NMR	Conform	Report attached

**Storage condition:** Store at 2-8°C, protect from direct sunlight.**Shipment date:**

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

**Manufacturing date:****Retest date:****Document data reference:***Prepared by**Reviewed by**Approved by*