

CERTIFICATE OF ANALYSIS


Reference Substance

Date

Dipyridamole Impurity E

Batch No.

IUPAC Name: 2,2',2'',2'''-[[6,8-di(piperidin-1-yl)pyrimido[5,4-d]pyrimidine-2,4-diyl]dinitrilo]tetraethanol

Molecular Formula: C₂₄H₄₀N₈O₄

Molecular weight: 504.63g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Yellow Colour solid	-----
2	Purity by HPLC	98.04%	Chromatogram attached
3	Mass	Conform	Report attached
4	NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight.**Shipment date:**

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest date:****Document data reference:***Prepared by**Reviewed by**Approved by*