

CERTIFICATE OF ANALYSIS


Reference Substance

Date

PPG ESTER-2

Batch No.:

IUPAC Name: [2-(2, 6-Dichloro-phenylamino)-phenyl]-acetic acid 2-hydroxy-1-methyl-ethyl ester

Molecular Formula: C₁₇H₁₇Cl₂NO₃

Molecular weight: 354.2

Sr.#	Test	Results	Data ref.
1	Appearance	Off-white solid	-----
2	Purity by HPLC	74.02%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store in refrigerator and protect from moisture**Shipment date:**

This certificate of analysis is valid for one year from date of shipment, provided the substance is stored under suitable conditions.

Retest date:**Document data reference:***Prepared by**Reviewed by**Approved by*