

CERTIFICATE OF ANALYSIS


Reference Substance

Date:

Diclofenac Impurity A

Batch No.

IUPAC Name: 1-(2,6-dichlorophenyl)-1,3-dihydro-2H-indol-2-one.

Molecular Formula: C₁₄H₉Cl₂NO

Molecular weight: 278.13

Sr.#	Test	Results	Data ref.
1	Appearance	Off white Solid	-----
2	Purity by HPLC	99.50%	Chromatogram attached
3	¹ H NMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store at 2-8°C, Protect from direct sunlight. Protect from moisture.**Shipment date:**

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest date:****Document data reference:***Prepared by**Reviewed by**Approved by*