

CERTIFICATE OF ANALYSIS

Reference Substance

Date:


Diclofenac Impurity F

Batch No.

IUPAC Name: N-Chloroacetyl-(2,6-dichlorophenyl)phenylamine

OR

(Alpha-(2,6-Dichlorophenyl)-4-chloroacetanilide)

Molecular Formula: C₁₄H₁₀Cl₃NO

Molecular weight: 314.0

Sr.#	Test	Results	Data ref.
1	Appearance	A light purple colour Solid	-----
2	Purity by HPLC	99.06%	Chromatogram attached
3	¹ H NMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store at 2-8°C, Protect from direct sunlight. Protect from moisture.**Shipment date:**

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest date:****Document data reference:***Prepared by**Reviewed by**Approved by*