

## CERTIFICATE OF ANALYSIS


Reference Substance

Date:

**1-(2,6-Dichlorophenyl)isatin**

Batch No.:

IUPAC Name: 1-(2,6-Dichlorophenyl)isatin

Molecular Formula: C<sub>14</sub>H<sub>7</sub>Cl<sub>2</sub>NO<sub>2</sub>

Molecular Weight: 292.00 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Yellow colour Solid	-----
2	Purity by HPLC	99.54%	Chromatogram attached
3	MASS	Conform	Report attached
4	<sup>1</sup> H NMR	Conform	Report attached

**Storage condition:** Store at 2-8 °C, protect from direct sunlight and Moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

**Manufacturing date:****Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*