

## CERTIFICATE OF ANALYSIS


Reference Substance

Date:

**2-(2,6-Dichloroanilino)phenylglyoxylic acid/Diclofenac glyoxilic acid**

Batch No.:

IUPAC Name: 2-(2,6-Dichloroanilino)phenylglyoxylic acid

Molecular Formula: C<sub>14</sub>H<sub>8</sub>Cl<sub>2</sub>NO<sub>3</sub>Na

Molecular Weight: 308.00 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Yellow colour Solid	-----
2	Purity by HPLC	99.70%	Chromatogram attached
3	MASS	Conform	Report attached
4	<sup>1</sup> H NMR	Conform	Report attached

**Storage condition:** Store at 2-8 °C, protect from direct sunlight and Moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

**Manufacturing date:****Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*