

CERTIFICATE OF ANALYSIS


Reference Substance

Date-

Nebivolol Related Compound-B

Batch No.:

IUPAC Name: [2S*[1R*, 5S*(S*)]]-alfa,alfa'-[Iminobis(methylene)]bis[6-fluoro-3,4-dihydro-2H-1-benzopyran-2-methanol]


 Molecular Formula: C₂₂H₂₅F₂NO₄

Molecular weight: 405.44 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	White Solid	-----
2	Purity by HPLC	99.82%	Chromatogram attached
3	Mass	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Storage condition: Store at 2-8°C, Protect from direct sunlight and moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:
Retest date:
Document data reference:
Prepared by
Reviewed by
Approved by