


CERTIFICATE OF ANALYSIS

Reference Substance

Date:

Nebivolol Desfluoro (Free base)**Batch No.:**

IUPAC Name: 1-(chroman-2-yl)-2-(2-(6-fluorochroman-2-yl)-2-hydroxyethylamino)ethanol

Molecular Formula: C₂₂H₂₆FNO₄

Molecular Weight: 387.4 g/ mol

	Test	Results	Data ref.
1	Appearance	Off White solid	-----
2	Purity by HPLC	99.14%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached

Storage condition: Store at 2-8 degree, protect from direct sunlight and moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*