

CERTIFICATE OF ANALYSIS


Reference Substance

Date:

Montelukast other-Isomer Impurity

Batch No.:

IUPAC Name: (S,E)-2-(1-(((1-(3-(2-(7-chloroquinolin-2-yl)vinyl)phenyl)-3-(2-(2-hydroxypropan-2-yl)phenyl)propyl)thio)methyl)cyclopropyl)acetic acid


 Molecular Formula: C₃₅H₃₆ClNO₃S

Molecular weight: 586.18g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off White Solid powder	-----
2	Purity by HPLC	99.48%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:

Retest date:

Document data reference:

Prepared by

Reviewed by

Approved by