

CERTIFICATE OF ANALYSIS

Reference Substance


Date:

Acetylsalicylsalicylic Acid / Aspirin EP Impurity D**Batch No.:**

IUPAC Name: Acetylsalicylsalicylic Acid

OR

(2-[2-acetylhydroxybenzoyl] oxy) benzoic Acid)

Molecular Formula: $C_{16}H_{12}O_6$

Molecular Weight: 300.26 g/mol.

Sr.#	Test	Results	Data ref.
1	Appearance	White colour solid	-----
2	Purity by HPLC	99.63%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached

Storage condition: Store at 2-8 degree, protect from direct sunlight and moisture.**Shipment Date:**

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*