

## CERTIFICATE OF ANALYSIS


Reference Substance

Date:

**Naproxen Impurity H**

Batch No.:

IUPAC Name: 6-methoxynaphthalen-2-ol



Molecular Weight: 174.20 g/mol

Molecular Formula: C<sub>11</sub>H<sub>10</sub>O<sub>2</sub>

| Sr.# | Test               | Results           | Data ref.             |
|------|--------------------|-------------------|-----------------------|
| 1    | Appearance         | Light Brown Solid | -----                 |
| 2    | Purity by HPLC     | 96.43%            | Chromatogram attached |
| 3    | Mass               | Conform           | Report attached       |
| 4    | <sup>1</sup> H NMR | Conform           | Report attached       |

**Storage condition:** Store at 2-8°C, protect from moisture and direct sunlight.

This certificate of analysis is valid for two years from date of manufacturing provided the substance is stored under suitable conditions.

**Date of manufacturing:****Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*