

CERTIFICATE OF ANALYSIS

Reference Substance


Date:

NORETHINDRONE ACETATE IMP-C (EP)

Batch No.:

IUPAC Name: Δ -5(6)-Norethindrone Acetate;

OR

17 β -Acetoxy-19-nor-17 α -pregn-5-en-20-yn-3-oneMolecular Formula: $C_{22}H_{28}O_3$

Molecular weight: 340.46g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	White solid	-----
2	Purity by HPLC	82.99%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store at 2-8°C, Protect from direct sunlight and moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest date:****Document data reference:***Prepared by**Reviewed by**Approved by*