


## CERTIFICATE OF ANALYSIS

Reference Substance

Date:

**Fluconazole Related Compound B****Batch No.:**

IUPAC Name: (2-(4-fluorophenyl)-1,3-bis(1H-1,2,4-triazol-1-yl)-propan-2-ol)

Molecular Formula: C<sub>13</sub>H<sub>13</sub>FN<sub>6</sub>O

Molecular Weight: 288

Sr.#	Test	Results	Data ref.
1	Appearance	A off White Colour Solid	-----
2	Purity by HPLC	92%	Chromatogram attached
3	<sup>1</sup> HNMR	Conform	Report attached
4	MASS	Conform	Report attached

**Storage condition:** Store at 2-8 degree, protect from direct sunlight.**Shipment date:**

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

**Manufacturing date:****Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*