


CERTIFICATE OF ANALYSIS

Reference Substance

Furosemide Impurity-D

Batch No.:

IUPAC Name: 2, 4-bis [(furan-2-ylmethyl) amino]-5-sulphamoylbenzoic acid.



Molecular Formula = C₁₁H₁₁N₃O₆S
Formula Weight = 391.39838

Sr.#	Test	Results	Data ref.
1	Appearance	White to off white	-----
2	Purity by HPLC	>95%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store at room temperature in desiccator and protect from moisture.

This certificate of analysis is valid for one year from date of shipment provided the substance is stored under suitable conditions.

Document data reference:

Prepared by

Reviewed by

Approved by