

CERTIFICATE OF ANALYSIS


Reference Substance

Date:

Trazodone Impurity B

Batch No. :

IUPAC Name: 2-(3-(4-phenylpiperazin-1-yl)propyl)-[1,2,4]triazolo[4,3-a]pyridin-3(2H)-one

Molecular Formula: C₁₉H₂₃N₅O

Molecular weight: 337.42 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off White Solid	-----
2	Purity by HPLC	95.48%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest date:**

Document data reference:

Prepared by

Reviewed by

Approved by