

CERTIFICATE OF ANALYSIS


Reference Substance

Date:

Trazodone Impurity -G

Batch No.

IUPAC Name: 1-(3-chlorophenyl)-4-(3-isobutoxypropyl)piperazine hydrochloride

Molecular Formula: C₁₇H₂₈Cl₂N₂O

Molecular weight: 347.32g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Light Brown solid	-----
2	Purity by HPLC	86.22%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight and Moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest date:**

Document data reference:

*Prepared by**Reviewed by**Approved by*