

## CERTIFICATE OF ANALYSIS


Reference Substance

Date:

**Gliclazide impurity F**

Batch No.

IUPAC Name: N-((hexahydrocyclopenta[c]pyrrol-2(1H)-yl)carbamoyl)-2-methylbenzenesulfonamide

Molecular Formula: C<sub>15</sub>H<sub>21</sub>N<sub>3</sub>O<sub>3</sub>S

Molecular weight: 323.41

Sr.#	Test	Results	Data ref.
1	Appearance	White solid	-----
2	Purity by HPLC	96.47%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform	Report attached

**Storage condition:** Store at 2-8°C, protect from direct sunlight and Moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

**Manufacturing date:****Retest date:****Document data reference:***Prepared by**Reviewed by**Approved by*