


CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard

Date:

Gliclazide EP impurity B**Catalogue no.: VL650002****Batch No.:****IUPAC Name:** 2-Nitroso-octahydrocyclopenta[c]pyrrole; 3-Nitroso-3-azabicyclo[3.3.0]octaneMolecular Formula: C₇H₁₂N₂O

Molecular weight: 140.18 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	White solid	-----
2	Purity by HPLC	99.62%	Chromatogram attached
3	Mass	Conform	Report attached
4	NMR	Conform	Report attached

Material Shipping condition: At Room Temperature**Long Term Storage:** Store at 2-8 degree, Protect from Direct Sunlight and Moisture.

This certificate of analysis is valid for three years from date of manufacturing provided the substance is stored under suitable conditions.

Manufacturing date:**Retest date:****Document data reference:***Prepared by**Reviewed by**Approved by*