

CERTIFICATE OF ANALYSIS


Reference Substance

Date

Glipizide Related Compound A

Batch No.

IUPAC Name: 5-methyl-N-[2-(4-sulphamoylphenyl) ethyl] pyrazine-2-carboxamide

Molecular Formula: C₁₄H₁₆N₄O₃S

Molecular weight: 320.37g/mol

| Sr.# | Test | Results | Data ref. |
|------|----------------|-----------------|-----------------------|
| 1 | Appearance | Off White solid | ----- |
| 2 | Purity by HPLC | 99.5% | Chromatogram attached |
| 3 | MASS | Conform | Report attached |
| 4 | NMR | Conform | Report attached |

Storage condition: Store at 2-8 degree temperature**Shipment date:**

This certificate of analysis is valid for one year from date of shipping, provided the substance is stored under suitable conditions.

Retest date:**Document data reference:***Prepared by**Reviewed by**Approved by*