


CERTIFICATE OF ANALYSIS

Reference Substance
Domperidone Impurity A
Batch No.:

Date:

IUPAC Name: 5-chloro-1-(piperidin-4-yl)-1H-benzo[d]imidazol-2(3H)-one.



Molecular Weight: 251.71 g/mol
Molecular Formula: C₁₂H₁₄ClN₃O

Sr.#	Test	Results	Data ref.
1	Appearance	White Solid powder	-----
2	Purity by HPLC	99.90%	Chromatogram attached
3	MASS	Conform	Report attached
4	NMR	Conform	Report attached

Storage condition: Store at 2-8°C and protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by