


CERTIFICATE OF ANALYSIS

Reference Substance

Date:

Domperidone Impurity F**Batch No.:****IUPAC Name:** 1,3-bis[3-[4-(5-Chloro-2-oxo-2,3-dihydro-1H-benzimidazol-1-yl)piperidin-1-yl]propyl]-1,3-dihydro-2H-benzimidazol-2-one.

Molecular Weight: 717.79 g/mol

Molecular Formula: C₃₇H₄₂Cl₂N₈O₃

Sr.#	Test	Results	Data ref.
1	Appearance	Orange Colour Solid	-----
2	Purity by HPLC	97.59%	Chromatogram attached
3	MASS	Conform	Report attached
4	NMR	Conform	Report attached

Storage condition: Store at 2-8°C and protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*