

CERTIFICATE OF ANALYSIS


Reference Substance

Date:

Metoclopramide Impurity E/ N,N-diethylethane-1,2-diamine

Batch No.:

IUPAC Name: N,N-diethylethane-1,2-diamine.



Molecular Weight: 116.20 gm/ mol

Molecular Formula: C₆H₁₆N₂

Sr.#	Test	Results	Data ref.
1	Appearance	Pale yellow liquid	-----
2	Purity by GC	94.75%	Chromatogram attached
3	Mass	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight & Moisture

This certificate of analysis is valid for three years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*