

CERTIFICATE OF ANALYSIS


Reference Substance

Date:

Clarithromycin Impurity F

Batch No.:

IUPAC Name: 12-O-Methyl Clarithromycin; 6,12-Di-O-methylerythromycin



Molecular Weight: 761.98 g/mol

Molecular Formula: C₃₉H₇₁NO₁₃

Sr.#	Test	Results	Data ref.
1	Appearance	White solid	-----
2	Purity by HPLC	92.43%	Chromatogram attached
3	Mass	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight & Moisture

This certificate of analysis is valid for three years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by