

CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD


Reference Standard

Date:

Warfarin EP Impurity C/ Benzalacetone/ Benzylideneacetone**Catalog NO: VL910003**

Batch No.:

IUPAC Name: (3E)-4-Phenylbut-3-en-2-one;



Molecular Weight: 146.19 g/mol

Molecular Formula: C₁₀H₁₀O

| Sr.# | Test | Results | Data ref. |
|------|--------------------|------------------|-----------------------|
| 1 | Appearance | Pale Beige Solid | ----- |
| 2 | Purity by HPLC | >95.00% | Chromatogram attached |
| 3 | Mass | Conform | Report attached |
| 4 | ¹ H NMR | Conform | Report attached |

Storage condition: Store at 2-8°C, protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*