

## CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard


Date:

**Eletriptan Related Compound 01**

**Catalogue no.: VL990003**

**Batch No.:**

IUPAC Name: 3-[[[(2S)-1-Methyl-2-pyrrolidinyl]methyl]-5-[2-(phenylsulfonyl)ethyl]-1H-indole



Molecular Formula: C<sub>22</sub>H<sub>26</sub>N<sub>2</sub>O<sub>2</sub>S.HBr

Molecular Weight: 382.52 g/ mol (Free Base)

Sr.#	Test	Results	Data ref.
1	Appearance	Off White Solid	-----
2	Purity by HPLC	99.51%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached

**Storage condition:** Store at 2-8 degree, protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

**Manufacturing date:**

**Retest Date:**

**Document data reference:**

*Prepared by*

*Reviewed by*

*Approved by*