


CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard

Date:

Eletriptan Related Compound 02**Catalogue no.:** VL990002**Batch No.:**

IUPAC Name: 3-((R-1-Methylpyrrolidin-2-yl)methyl)-5-(E-2-(phenylsulfonyl)vinyl)-1H-indole

Molecular Formula: C₂₂H₂₄N₂O₂S

Molecular Weight: 380.5 g/ mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off white solid	-----
2	Purity by HPLC	99.29%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached

Storage condition: Store at 2-8 degree, protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*