


CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard

Date:

Eletriptan related compound 03**Catalogue no.: VL990001****Batch No.:**

IUPAC Name: 1H-indole,5-ethyl-3-[[[(2R)-1-methyl-2-pyrrolidinyl]methyl]

Molecular Formula: C₁₆H₂₂N₂

Molecular Weight: 242.36 g/ mol

Sr.#	Test	Results	Data ref.
1	Appearance	Dark Brown sticky	-----
2	Purity by HPLC	91.61%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached

Storage condition: Store at 2-8 degree, protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*